



DISPOSABLE
FORCEPS



VECTEC
SINGLE USE
ENDOSCOPIC DEVICES



VECTEC

VECTEC MONOPOLAR SINGLE USE LAPAROSCOPIC GRASPING FORCEPS



REF : * MD5 330
* MD5 470



REF : F5 330, FS5 330
FS5 470.



REF : ECS5 330, ECSS 470.

Reference	Type of jaws	Working length	Diameter	Jaws size
* MD5 330	● Dissector Maryland	330 mm	5 mm	20 mm
* D5 330	● Dissector	330 mm	5 mm	18 mm
A5 330	● Atraumatic	330 mm	5 mm	18 mm
ECS5 330	● Endoclinch	330 mm	5 mm	22 mm
F5 330 FS5 330	● Fenestrated	330 mm 330 mm	5 mm	18 mm
G5 330	● Claw, 2 x 6 teeth	330 mm	5 mm	18 mm

* without ratchet ● Yellow: Dissector ● Green: Atraumatic ● Purple: Traumatic

Ergonomic handle with ratchet mechanism guarantees strong and precise control with minimal effort during grasping and holding tissues.

Double-action jaws allows perfect grip.

Jaws rotates 360° and are adjusted using the rotational knob at the front of the handle.

The VECTEC grasping forceps comes fitted with a standard monopolar connector pin.

Electrical insulation is total right up to the active portion of the jaws.

The VECTEC grasping forceps is packed in a sealed pouch and sterilized by ETO to EN ISO 11135-1 standards.

Packaging : box of 5 individually packaged units.

MD Class IIa and IIb

Made in France

See instructions for use

Cleared for marketing by FDA



VECTEC

Bioparc - Vichy - Hauterive
03270 Hauterive
FRANCE

Phone : +33 (0)4 70 58 41 70
Fax : +33 (0)4 70 58 41 79
Email : vectec@vectec.fr
<http://www.vectec.fr>