

RECTAL CANNULA



**VECTEC**  
SINGLE USE  
ENDOSCOPIC  
DEVICES

**VECTEC**



## Single-use VECTEC cannula

- Allows excellent recto sigmoid exposure.
- Facilitates double safety testing.
  - 1- Air test
  - 2- Liquid dye test.



Reference	Désignation	Shaft diameter	Tip diameter	Tip length	No. unit per box
RCC 320	Rectal cannula	4 mm	32 mm	50 mm	10

### product description :

- A distal cannula with:
  - 1- an ergonomic and atraumatic tip for safe progression into the rectal and sigmoid cavities.
  - 2- a flat rear part to the tip good for safety tests (air/liquid).
- A small-diameter stainless steel shaft which is flexible to adapt to the patient's anatomy.
- An ergonomic handpiece.
- T tubing:
  - 1- a length of small diameter silicone tubing with a stop-cock valve and luerlock for air/liquid injection, to be located within the sterile zone.
  - 2- a length of large diameter tubing with a red clamp for waste.



VECTEC certifications :  
 EN ISO 9001 (2008)  
 EN ISO 13485 (2012)  
 ISO 13485 (2003) - CMDCAS

The VECTEC rectal cannula is packed in a sealed pouch and ETO sterilized to NF EN ISO 11135-1 standards.

MD Class IIa  
 See instructions for use  
 Made in France

CE  
 0197

© G R A  
 P H I C A  
 0 6 8 0 6  
 8 4 9 7 4